

M.A.S. Medical
For Your All of Your Veterinary Needs

D.B.I. Products

Helping to Improve Patient Dental Health

VE-2N

All-in-one Dental Unit

Part# VE-2N Delivery Unit with Suction
Part# VE-2NL Delivery Unit with Curing Light

5K Polisher

Piezo Ultrasonic Scaler with LED

Suction System

Swivel Lighted Push Button High Speed Hand piece with Quick Disconnect

Air Water-Spray Syringe

Two Quick Release Reservoirs (2 litres each)

Main Features:

- Automatic Instrument Selection
- Quiet Oil Free Compressor 6L, 1.2 HP
- Hair Catch Filters Increasing Compressor's Life Span
- 360 Swivel Fibre Optic Push Button High Speed Hand piece with Quick Disconnect Coupler
- Piezo Ultrasonic Scaler with 70,000 LUX of Light
- Scaler Accessories Include:
 - 3 Scaler Tips, 3 Sleeves, 3 Torque Wrenches
- Starter Kit Includes:
 - Bur Set, Suction Tips (HVE and Softflex with Adaptor), Hand piece Lubricant, Suction System Cleaningsolution, Disposable Prophylaxis Angles + Polishing Paste
- Quick and Hassle Free Set-Up
- 2 Year Limited Warranty

1.25 Horse Power oil Free Compressor
(specially designed for Veterinary Applications)

4 Medical Grade Caster Wheels with Brakes

Suction Collection Bottle

Integrated Foot Control (for all instruments)



M.A.S. Medical

For Your All of Your Veterinary Needs

SELF CONTAINED DENTAL UNIT

VE-4



Main Features:

- Automatic Instrument Selection
- Quiet Oil Free Compressor 4L, 3/4 HP
- 360 Swivel Lighted Push Button High Speed Hand piece with Quick Disconnect Coupler
- Piezo Ultrasonic Scaler without LED
- Scaler Accessories Include:
3 Scaler Tips, 3 Sleeves,
3 Torque Wrenches
- Starter Kit Includes:
Bur Set, Hand piece Lubricant,
Disposable Prophy Angles
- Quick and Hassle Free Set-Up
- 2 Year Limited Warranty

All-in-one Dental Unit

Part# VE-4

Air-Water-Spray Syringe

5k Polisher

Piezo Ultrasonic Scaler without LED

Swivel Lighted Push Button high Speed hand piece with Quick Disconnect

One Quick Release Reservoir

1.25 Horse Power oil Free Compressor (specially designed for Veterinary Applications)

4 Medical Grade Caster Wheels with Brakes

Suction Collection Bottle

Integrated Foot Control (for all instruments)

A Wide Variety of Products At Exceptional Prices

Handpieces



20E Air Motor



Straight Nose



5 SGII Low Speed Polisher



20 SGII Low Speed Handpiece



Contra Angle



Maxi Torque



Maxi Torque Fiber Optic



JP390



**Free Dental Stool
With the Purchase
the VE-4**

(retail \$475.50)

Buy Two, Get One Free! Limited Time Offer.



M.A.S. MEDICAL SUPPLIES INC.
MEDICAL • DENTAL • VETERINARY

Contact us for more information:

36 Lauraview Crescent, Unionville, Ontario, L6C 2R4

www.masmedicalsupplies.com

Email: sales@masmedicalsupplies.com Phone: (905) 927-0548 Fax: (905) 887-5306

VISA

