

M.A.S. Medical

For All Your Veterinary Needs
Helping Improve Patient Dental Health

EFFECTIVE • ROBUST • PRACTICAL



VE-4M

- Compact, Strong and Reliable
- Premium Hand-piece Package
- Piezo Scaler w/LED Illumination
- All Enclosed, 3/4 HP, Oil-Free Compressor
- Integrated Foot Control for All Instruments

(Full specs. on reverse)

8 Year Warranty Protection Plan
– Best in the Business!



VE-2ME

- Premium Hand-piece Package
- Piezo Scaler w/LED Illumination
- Quiet, 1.2 HP, Oil-Free Compressor
- AVS Suction System
- Integrated Foot Control for All Instruments

(Full specs. on reverse)

8 Year Warranty Protection Plan
– Best in the Business!

TT-6 Plus Stand

- Compact, Lightweight, Self-Contained
- Stainless Steel Stand w/ a 24" x 24" 5-Caster Base
- Height Adjustable up to 40"
- All Enclosed, 1/3 HP, Oil-Free Compressor
- Integrated Foot Control for All Instruments

(Full specs. on reverse)

8 Year Warranty Protection Plan
– Best in the Business!



**FREE DENTAL STOOL
WITH PURCHASE**

(retail value \$475.50)



FEATURES	VE-2ME	VE-4M	TT-6
HIGH & LOW SPEED HANDPIECES			
Swivel High Speed, Push Button, Fibre Optic Handpiece (Drill) with up to 400,000 rpm	•	•	OPT
Fixed High Speed, Push Button, Fibre Optic Handpiece (Drill) with up to 350,000 rpm	N/A	N/A	•
Polisher Low Speed Handpiece; Reduction 4:1 (More Torque), 5,000 rpm; Forward/Reverse Direction	•	•	•
PIEZO SCALER			
Piezo Ultrasonic Scaler with a Detachable, Aoclavable Handpeice with Direction Adjustable LED	•	N/A	N/A
Piezo Ultrasonic Scaler, Detachable, all Metal Handpiece with LED	N/A	•	•
Three (3) Piezo Ultrasonic Scaler Tips Included	•	•	•
WATER CAPACITY			
Dual Water Reservoir of 2 Liter Capacity each w/ Quick Disconnect	•	N/A	N/A
Singular Water Reservoir of 2 Liter capacity w/ Quick Disconnect	N/A	•	N/A
Singular Water Reservoir of 750ml Capacity	N/A	N/A	•
COMPRESSORS			
Enclosed, Oil-Free Air Compressor at 1.2 HP with Hair Catch Filter plus Regular Internal Filter, Automatic Electric Cooling Fan, Overload Protection Switch, Noise Level at 57db, Sound Insulated	•	N/A	N/A
All Enclosed, Oil-Free Air Compressor at 3/4 HP with a Regular Internal Filter, Automatic Electric Cooling Fan, Noise Level at 55db	N/A	•	N/A
All Enclosed, Oil-Free Air Compressor at 1/3 HP with a Regular Internal Filter, Automatic Electric Cooling Fan, Noise Level at 53db	N/A	N/A	•
Compressor designed with a Hassle Free Purge Valve with Easy Access	•	•	•
Fuse Protection for Compressor	•	•	•
Compressor equipped with Internal Safety Valve	•	•	•
Nitrogen Connection in lieu of Air Compressor	OPT	OPT	OPT
SPECIFICATIONS			
Maximum Height (from Top Tray)	37"	37.5"	43"
Minimum Height (from Top of Tray)	27.75"	26.5"	31.75"
Weight Including Compressor (approx.)	104 lbs	90 lbs	85 lbs
Base Frame Width	18"	18"	23"
Dental Unit Head Size (Length, Width, Height)	L: 20" W: 14.5" H: 4.75"	L: 14" W: 14" H: 3.5"	L: 12" W: 8" H: 4.5"
ADDITIONAL FEATURES			
Suction System with Collection Bottle activated at Handpiece Lever	•	N/A	N/A
Master Valve to Disconnect Air and Water. Refill Reservoir Safely, Protects unit after hour	•	•	•
Additional Water Port for other Equipment additions	•	•	N/A
Water On/Off Switch	•	•	•
Instrument Pressure Gauge	•	•	•
Instrument Metal Tray	•	•	OPT
6 Outlet Power Strip w/ Surge Protection, with On/Off Switch	•	•	•
Curing Light, 1,500mw/cm ² (OPT)	OPT	OPT	OPT
Integrated Foot Control that Activates all Instruments, Metal Surface, Circular Form for Easy Activation	•	•	•
Three Way Syringe w/ Quick Release Autoclavable Metal Tip	•	•	•
Starter Kit includes: (10 Assorted Burs, Lubricant, Disposable Prophy Angles w/ Sample Paste, (*) High Volume Ejector and SoftFlex Suction Tips, 1 Adaptor, Suction System Cleaning Solution.	•	•	•

(*) For Dental Units with Suction Systems | TBD = To Be Determined by Client's choice



Contact us for more information:

36 Lauraview Crescent, Unionville, Ontario, L6C 2R4

www.masmedicalsupplies.com

Email: sales@masmedicalsupplies.com Phone: (905) 927-0548 Fax: (905) 887-5306

