

Trust M.A.S. Medical for all your Surgical Equipment, Stainless Steel Cages and Accessories

Medical Lighting Surgical

Model LED 130F Focusable



The VSSI LED-130F focusing function is activated by turning the handle on the light head.

- Focusable light field: 5.5" - 9.8"
- Light intensity at a distance of 39.4": 70,000 lux or 6,503 foot-candles
- Color temperature: 4,300 kelvin
- Light head diameter: 13"
- Number of LEDs: 19
- Anticipated life of LEDs: >40,000 hours
- 3-year limited warranty
- sterilizable handle available

Mounting Options

LED-120 or LED-130 Series Lights



Medical Lighting Exam Light

Model LED 120

- Fixed-focus light field: 5.1"
- Light intensity at a distance of 39.4": 35,000 lux or 3,252 foot-candles
- Color temperature: 4,300 kelvin
- Light head diameter: 11.4"
- Number of LEDs: 12
- Anticipated life of LEDs: >40,000 hours
- 3-year limited warranty
- sterilizable handle available



Stainless Steel Cages

Easy To Assemble
Completely Interchangeable

Easier, More Stable Assembly

VSSI's unique "Lap-Lock" overlapping, locking, uniframe mounting system makes assembly much easier and the overall structure more durable.

Multiple 6" centered self-tapping holes within the uniframe system are utilized to securely lock cages together.



Call For
Package
Pricing

Shown With Accessories:
Side Panel, Top Panel
and Mobile Platform



Call for Custom Pricing!

Cost Varies By Number of Units, Size, and Configuration

Plastisol Coated Cage Floors

Designed to give added
comfort to patients



They are designed to give added comfort to patient and will keep the animal warm, dry and clean.


Available in 18" / 24" / 30" / 36" and 48" cage or kennel floors. *

These unique Raised floors will fit all of VSSI / Shore Line and Suburban models.

***Note:** 48" Cages with centre divider require two pieces

Some Popular Cage Configurations

Many other configurations available.



- A Top Row: Two 24" x 30"
Bottom Row: One 48" x 30" Double Door
- B Three 24" x 30"
Bottom Row: Two 36" x 30"
- C Top Row: Four 30" x 24"
Bottom Row: One 48" x 30" Double Door, Two 36" x 30"

- D Top Row: Four 24" x 24"
Middle Row: Two 36" x 30", One 24" x 30"
Bottom Row: Two 48" x 30"
- E Top Row: Five 24" x 24"
Middle Row: Two 36" x 30", One 24" x 30"
Bottom Row: One 48" x 30" Double Door, Two 36" x 30"

Peninsula™ Folding Lift Table

Lifts ■ Weighs ■ Folds
(optional feature)



Restraint System

Three-Sided Work Area

22" x 44"
Stainless Steel Top With
Raised Anti-Drip Edges

Shown With Optional Built-In
The VSSI Way™ Scale

Depress Handle
To Fold Table

225 Pound
Lift Capacity

Height Adjusts
8 1/2" To 40"
From Floor
(scale top model)



Choose Hand or Foot Control

Folds At
Any Height
To 16" From Wall



Peninsula Lift, Hand Control
Peninsula Lift, Foot Control
Peninsula Lift With Scale, Foot Control

Premier Mobile Lift Table

22" x 50"
Stainless Steel
Exam Surface With
Raised Anti-Drip Edges

Lifts Up To
225 Pounds

Adjusts From
12 1/2" To 37 1/4"

Heavy Duty
Sturdy Design



Rechargeable
Battery
Operated

Available With
Optional Scale Top

Premier Mobile Lift Table
Premier Mobile Lift Table With Scale
Rail Kit

Kennel Runs



Patent Pending



Standard Door

- Reversible stainless steel doors – side by side units now may be opened in **opposite directions**.
- Stainless steel door hinges are mounted on top and bottom – not the side.

Sliding Door

- ★ **INDUSTRY FIRST** ★
- Convenient one-hand operation
- Maximizes kennel room(s) eliminates the need for "swing door" space
- Maintenance free marine grade anodized aluminum glide
- Also available with glass panels

Glass Door

- Open view boarding for owner comfort or full view observation
- Easy, convenient cleaning
- Aesthetically pleasing
- Scratch resistant
- Also available as sliding model



M.Q.S. MEDICAL SUPPLIES INC.
MEDICAL • DENTAL • VETERINARY

36 Lauraview Crescent, Unionville, Ontario, L6C 2R4

www.masmedicalsups.com

Email: sales@masmedicalsups.com Phone: (905) 927-0548 Fax: (905) 887-5306